

LT-R28A User's Manual

Concept:

The Effect of Light on Cell: Increase collagen synthesis and cellular proliferation

There are photoreceptors at the molecular - cellular level which, when triggered, activate a number of biological reactions: DNA / RNA synthesis, increase protein and collagen synthesis, and cellular proliferation. The result is rapid regeneration, normalization and healing of damaged cellular tissue. In essence, light is a trigger for the rearrangement of cellular metabolism.





A. Machine:

- Main Machine
- Stand

B. Parts & Accessories:

No.	Descriptions	Qty.
1	Colored Light Cover	1PCE





Color's Functions:

Color	Descriptions	Applications
Red	Stimulates us and always appears when the circulation is in full swing. Red, the warmest color, symbolizes life and fills us with strength and energy.	Vitalize Improve circulation Facilitate the flow of blood
Blue	Has a calming effect and helps us sleep better. Blue is the coldest, purest and deepest color and represents reserve and introversion.	Calm Increase synthesizing of protein and collagen Useful to sensitive, oily skin and pimples. Help to Slimming
Green	Stands for natural vigour and is the most calming color, symbolizing hope and satisfaction. Green is synonymous with spiritual balance and harmony.	Balance Normalize sebaceous glands Reduce black heads and sores
Purple	Restrain the desire to eat. Has the effect of draining lymphatic system. It is a combination of red and blue.	Drain Renew skin, reduce wrinkles
Orange	Makes us feel cheerful and generates the body, activating the cell tissue and metabolism. It is associated with extrovertism, keeping	Brighten Whiten, decompose pigment, and smooth scars.

	melancholia at bay and giving us spiritual strength.	
Turquoise	Has a calming effect. It is a combination of blue and green.	Defense Strengthen the ability to protect of skin itself. Calm down, release stress, and cure sunburn
White	Guide skin to absorb	Guide Help skin to absorb nutrients, and anti-aging

Maintenance:

- Clean the handle with clean cloth (also can use little alcohol to sterilize) after each use.

Caution:

- Do not operate on patients who wear cardiac pacemaker.
- Do not operate on patients suffering from epilepsy, cancers, sunburn, serious sore, having surgery recently, and pregnant.

Attention:

- Store machines in a dry place
- Do not drop the machine or put heavy things on it
- Check the voltage before usage
- After using or during cleaning, please pull out the power cord
- Do not use wet cloth to wipe the panel

	R-28P
POWER SUPPLY	AC110V/220 Volts
FREQUENCY	50 Hz/60Hz
POWER CONSUMPTION	65Watts
FUSES	
DIMENSIONS	79*42.5*39cm
WEIGHT	Approx.17kgs